CANI-EX CRYO PSE, Freezing Media for Canine Sperm

Operating instructions

CANI-EX CRYO PSE is a frozen solution and consists of TRIS, citric acid, fructose, glutamine, LDL fragments, PolySE and glycerol @ 6.0% (V/V) in purified H₂O.

The general guidelines about sperm freezing like for example centrifugation time and speed for enrichment, target sperm concentration, straws filling, freezing methods and protocols are considered as known.

1. Take the frozen media (CANI-EX CRYO and CANI-EX) out of the freezer 24 hours in advance of use and place it in the fridge or take it out of the freezer 2 hours before use and place it in water bath at 30°C until unfrozen and afterwards to room temperature. If no frozen CANI-EX is available, it should be prepared if a wash procedure of the ejaculate needs to be done.
2. If the ejaculate is contaminated with prostatic fluid, a washing procedure is recommended. Measure or weigh the volume of the entire ejaculate. Dilute the ejaculate with the CANI-EX extender (1:3 V/V): meaning one part of semen and 3 parts of CANI-EX. Centrifuge for 10-15 minutes at 700 g. Remove the supernatant with a pipet and bring the volume back to the original volume of the ejaculate with CANI-EX Extender.
3. Measure the concentration, total motility, progressive motility and if needed the morphological parameters. Calculate the amount of CANI-EX CRYO extender that needs to be added.
4. Add slowly the CANI-EX CRYO extender to the 2nd fraction of the ejaculate or to the washed/enriched semen (*as described above*) until the desired concentration will be obtained.
5. Fill and seal the straws with the prepared diluted semen and keep the filled straws at +4°C for at least two hours.
6. Start the freezing procedure by a manual or by using a pre-programmed automatic freezer.
7. Transfer the frozen straws into a cannister into a liquid nitrogen container.

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